

Adapters an overview

23.2 MM PHOTO PORT C-MOUNT ADAPTERS

Suitable for all Ø 23.2 mm (photo) tubes

- EDC-1324** C-mount adapter with 0.45x objective for 1/2" camera
- EDC-1353** C-mount adapter with 0.37x objective, fixed focal length and short barrel*. For 1/3" camera
- EDC-1355** C-mount adapter with 0.50x objective, fixed focal length and short barrel*. For 1/2" camera
- EDC-1357** C-mount adapter with 0.75x objective, fixed focal length and short barrel*. For 2/3" camera
- EDC-1359** C-mount adapter with 1x objective, fixed focal length and short barrel, for standard 23.2 mm tube and 1 inch camera (no image available)

EDC-1330 30 mm to 23.2 mm conversion adapters

EDC-1335 30.5 mm to 23.2 mm conversion adapters



*for increased light efficiency

MODELS	1/3 inch camera	1/2 inch camera	2/3 inch camera	1 inch camera
EDC-1324		0.45x objective		
EDC-1353	0.37x objective, short barrel			
EDC-1355		0.5x objective, short barrel		
EDC-1357			0.75x objective, short barrel	
EDC-1359				1x objective, short barrel

Suitable for BioBlue, BioBlue.Lab, bScope, EcoBlue, iScope, MicroBlue and ZE-series

Delph-X Observer *, NexiusZoom * and StereoBlue *

Oxion * and Oxion Inverso * (all Oxion versions need extra EAE-5120-2 23.2 adapter)

DZ-series (needs extra DZ.9005 photo port attachment)

AP-series**, BE-series**, EduBlue** and P-series **

* For these models dedicated photo port adapters are also available

** These binocular microscopes need an extra EDC-1335 adapter

DIGITAL SLR ADAPTERS

Digital single lens reflex camera with a APS-C sensor can be mounted on a standard 23.2 mm trinocular tube. Choose the universal EDC-5130 adapter with built-in 2x objective, together with a brand-specific T-2 adapter

- EAE-5130** Universal SLR camera adapter for 23.2 mm tube with built-in 2x magnification lens. Needs T2 adapter
- EAE-5025** T2 adapter for Nikon D digital SLR cameras
- EAE-5040** T2 adapter for Canon EOS digital SLR cameras

Other T2 adapters on request



PHOTO PORT ADAPTERS

Photo port adapters work in a similar way to c-mount adapters but are shorter in design and of higher quality, resulting in more stable higher quality imaging.

Suitable for Euromex microscopes with a photo port such as Oxion, NexiusZoom, StereoBlue, Delphi-X Observer and Delphi-X Inverso



● EAE-5120-2

EAE-5120-2 Standard 23.2 mm diameter tube for Oxion photo port

EDI-9840 C-mount with high resolution 0.40x objective for 1/3 inch C-mount camera (no image available)

EDI-9850 C-mount with high resolution 0.50x objective for 1/2 inch C-mount camera (no image available)

EDI-9810 C-mount with 1x magnification for C-mount camera (no image available)

EDX-9835 C-mount with high resolution relay 0.35x objective for 1/3 inch C-mount camera

EDX-9850 C-mount with high resolution relay 0.50x objective for 1/2 inch C-mount camera

EDX-9863 C-mount with high resolution relay 0.63x objective for 2/3 inch C-mount camera

EDX-9810 C-mount with 1 magnification for C-mount camera

ENZ-9850 C-mount adapter with 0.5x objective for NexiusZoom microscopes and 1/2" cameras

EOX-9810 C-mount adapter with 1x magnification for 1" cameras

EOX-9833 C-mount adapter with 0.33x objective for Oxion microscopes and 1/3" cameras

EOX-9850 C-mount adapter with 0.5x objective for Oxion and Oxion Inverso microscopes and 1/2" cameras

ESB-9850 C-mount adapter with 0.5x objective for StereoBlue microscopes and 1/2" cameras



● EDX-9835

● EDX-9850

● EDX-9863

● EDX-9810



● ENZ-9850

● EOX-9810

● EOX-9833

● EOX-9850



● ESB-9850

MODELS	1/3 inch camera	1/2 inch camera	2/3 inch camera	1 inch camera
StereoBlue		ESB-9850		
NexiusZoom		ENZ-9850		
Oxion	EOX-9833	EOX-9850		EOX-9810
Delphi-X Observer	EDX-9835	EDX-9850	EDX-9863	EDX-9810
Delphi-X Inverso	EDI-9840	EDI-9850		EDI-9810

ADAPTERS FOR DIFFERENT BRANDS OF MICROSCOPES

MODELS	1/3 inch camera	1/2 inch camera	2/3 inch camera	1 inch camera
Leica ⁽¹⁾	EAE-9835-L (0.35x)	EAE-9855-L (0.55x)	EAE-9870-L (0.70x)	EAE-9810-L (1x)
Olympus ⁽²⁾	EAE-9835-O (0.35x)	EAE-9850-O (0.50x)	EAE-9863-O (0.63x)	EAE-9810-O (1x)
Zeiss ⁽³⁾	EAE-9835-Z (0.35x)	EAE-9850-Z (0.50x)	EAE-9865-Z (0.65x) EAE-9880-Z (0.80x)	EAE-9810-Z (1x)
Zeiss PrimoStar 3 ⁽⁴⁾	EAE-9835-ZP3 (0.35x)	EAE-9850-ZP3 (0.50x)	EAE-9865-ZP3 (0.65x)	EAE-9810-ZP3 (1x)
Zeiss ZA ⁽⁵⁾	--	EAE-9850-ZA (0.50x)	EAE-9863-ZA (0.63x) EAE-9880-ZA (0.80x)	EAE-9810-ZA (1x) EAE-9812-ZA (1.2x)
Nikon ⁽⁶⁾	EAE-9835-N (0.35x)	EAE-9855-N (0.55x)	EAE-9870-N (0.70x) EAE-9880-N (0.80x)	EAE-9810-N (1x)

⁽¹⁾ Leica DM, 35 mm to 25.4 mm C-mount, stainless steel

⁽²⁾ Olympus CX, CX2, BX, BX2, IX, MX, GX, SZX, STM, 42 mm to 25.4 mm C-mount, anodized aluminum

⁽³⁾ Zeiss PrimoStar, PrimoStar Vert, 30 mm to 25.4 mm C-mount, anodized aluminum

⁽⁴⁾ Zeiss PrimoStar3, C-mount, anodized aluminum

⁽⁵⁾ Zeiss 60N-C, 60N-T2 photo ports of Axio examiner, Axio Imager, Axio Lab, Axio Observer, Axio Scope, Axio Zoom, Stemy 508 and SteREO Discovery

⁽⁶⁾ Nikon E100, E200, 50i, 55i, 80i, 90i, Ni series, Ti Series, SMZ800, SMZ1000, SMZ15000, 38 mm to 25.4 mm C-mount, stainless steel

Standard C-Mount: 25.4 mm (1 inch) diameter with 32 threads per inch

LEICA

- For Leica microscopes with standard 35 mm outer diameter adapter (HC) to standard C-mount

EAE-9835-L 0.35x C-mount adapter for 1/3 inch cameras

EAE-9855-L 0.55x C-mount adapter for 1/2 inch cameras

EAE-9870-L 0.70x C-mount adapter for 2/3 inch cameras

EAE-9810-L 1x C-mount adapter for 1 inch cameras



OLYMPUS

- For Olympus microscopes with standard 35 mm outer diameter adapter to standard C-mount

- For Olympus models CX, CX2, BX, BX2, IX, MX, GX, SZX, STM


EAE-9835-O 0.35x C-mount adapter for 1/3 inch cameras

EAE-9850-O 0.50x C-mount adapter for 1/2 inch cameras


EAE-9863-O 0.63x C-mount adapter for 2/3 inch cameras

EAE-9880-O 0.80 C-mount adapter for 1 inch cameras


EAE-9810-O 1x C-mount adapter for 1 inch cameras




● EAE-9835-L




● EAE-9870-L




● EAE-9855-L




● EAE-9835-O



● EAE-9810-O



● EAE-9850-O





● EAE-9880-O

ZEISS P95 ADAPTER

- For Zeiss trinocular microscopes with standard Zeiss P95 adapter (30mm, 1.18 inch inner diameter phototube) to standard C-Mount
- For Zeiss models with UIS photo adapter such as PrimoStar (except PrimoStar 3), Primo vert

- EAE-9835-Z** 0.35x C-mount adapter for 1/3 inch cameras
EAE-9850-Z 0.50 C-mount adapter for 1/2 inch cameras
EAE-9865-Z 0.65x C-mount adapter for 2/3 inch cameras
EAE-9880-Z 0.80x C-mount adapter for 1 inch cameras
EAE-9810-Z 1x C-mount adapter for 1 inch cameras



ZEISS P90-C ADAPTER

- For Zeiss trinocular PrimoStar 3 microscopes with standard Zeiss P90 adapter

- EAE-9835-ZP3** 0.35x C-mount adapter for 1/3 inch cameras
EAE-9850-ZP3 0.50x C-mount adapter for 1/2 inch cameras
EAE-9865-ZP3 0.65x C-mount adapter for 2/3 inch cameras
EAE-9810-ZP3 1x C-mount adapter for 1 inch cameras



ZEISS 60N-T2

- For Zeiss trinocular microscopes with standard Zeiss 60N-T2 photo ports (60mm) to standard C-Mount
- For Zeiss models with 60N-T2 photo adapter such as Axio Examiner, Axio Imager, Axio Lab, Axio Observer, Axio Scope, Axio Vert.A1, Axio Zoom, Stemi 508 doc, Stemi 508 trino, SteREO Discovery, SteREO Lumar

EAE-9850-ZA 0.50x C-mount adapter for 1/2 inch cameras

EAE-9863-ZA 0.63x C-mount adapter for 2/3 inch cameras

EAE-9880-ZA 0.80x C-mount adapter for 2/3 inch cameras

EAE-9810-ZA 1x C-mount adapter for 1 inch cameras

EAE-9812-ZA 1.2x C-mount adapter for 1 inch cameras



NIKON

- For Olympus microscopes with standard 38 mm outer diameter adapter to standard C-mount
- For Nikon E100, E200, 50i, 55i, 80i, 90i, Ni series, Ti Series, SMZ800, SMZ1000, SMZ15000


EAE-9835-N 0.35x C-mount adapter for 1/3 inch cameras

EAE-9855-N 0.55x C-mount adapter for 1/2 inch cameras

EAE-9870-N 0.70x C-mount adapter for 2/3 inch cameras

EAE-9880-N 0.80x C-mount adapter for 2/3 inch cameras

EAE-9810-N 1x C-mount adapter for 1 inch cameras



www.euromex.com