


Digital stand-alone HDS

Euromex offers complete stand-alone camera systems with HD cameras which are suitable for each Euromex microscope. The solution contains an 11.6" LCD screen (which is directly connected to an HD camera) and offers you the option to take pictures, make videos, live view, measuring and lots of others options without the use of a personal computer! Just attach the USB mouse to your HD camera and you are ready to go!



VC.3020-HDS

HD-Mini camera, color HD high definition 1.2 MP high speed camera with built-on 11.6" HD screen, 1280 x 720p, stand-alone usage with standard 16 GB SD card, HDMI output. With built-in mouse-driven capture software. Can only be used in stand alone modus



VC.3036-HDS

HD-Ultra camera, color HD High definition high speed camera 1080p with built-on HD screen. 36 bit color rendering, 1/2.8 inch 6 MP Sony sensor, standard 16 GB SD card, HDMI and USB-2 output with built-in mouse-driven capture software. Can be used in stand alone or with a computer and supplied with ImageFocus Plus software. Supplied with 0.37x objective



VC.3031-HDS

HD-Lite camera, color HD High definition 5 MP high speed camera 1080p with HD screen. 24 bit color rendering, standard 8 GB SD card, HDMI and USB-2 output with With built-in mouse-driven software. Can be used in stand alone or with a computer and supplied with ImageFocus Plus software. Supplied with 0.45x objective



VC.3038-HDS

HD-Pro camera, color HD high definition 2.0 MP high speed camera 1080p with built-on 11.6" HD screen. 24 bits color rendering, 1/2.8 inch Sony sensor, standard 16 GB SD card, HDMI and USB-2 output with built-in mouse-driven capture software. Can be used in stand-alone mode or with a computer and is supplied with ImageFocus Alpha software. Comes with 0.37x objective



VC.3034-HDS

HD-Autofocus, color HD High definition high speed camera 1080p, 16 bit grey rendering, 36 bit color rendering, 1/1.9 inch Sony sensor, standard 16 GB SD card, HDMI output with built-in mouse-driven capture software. Can only be used in stand alone modus. Supplied with 0.50x objective



VC.3040-HDS

UHD-4K camera, color HD high definition 8.0 MP high speed camera 4K ultra HD with built-on 11.6" HD screen. 36 bits color rendering, 1/1.8 inch Sony 4K sensor, standard 16 GB SD card, HDMI, Ethernet, wifi and USB-3 output with built-in mouse-driven capture software. Can be used in stand-alone mode or with a computer or tablet and is supplied with ImageFocus Alpha software. Comes with 0.5x objective