

JuniorLED

HIGHLIGHTS

- High value for money
- Magnifications 40/100/400x
- Cordless LED illumination with internal batteries

TECHNICAL SPECIFICATIONS

EYEPIECES

Wide field WF 10x eyepieces

HEAD

Monocular 45° inclination, 360° rotatable

NOSEPIECE

For three objectives

OBJECTIVE

4x N.A. 0.10, 10x N.A. 0.25 objectives and a 540x N.A. 0.65 spring-loaded objective. Magnification 40x, 100x and 400x

STAND

Coarse adjustment with friction control

RACK STOP

Adjustable rack stop to prevent damage to sample and objectives.

STAGE

90 x 90 mm with two object clamps

CONDENSER

Single lens condenser with 6-positions diaphragm

ILLUMINATION

LED illumination with rechargeable batteries

PACKAGE CONTENT

In styrofoam with dust cover



● 81.500 JuniorLED

81.500 JUNIORLED

300 (h) x 120 (w) x 160 (d) | 1.6 kg