

Cold light source

Compact dual 2 W LED

HIGHLIGHTS

- Economical and compact light source
- For Materials Science en stereo microscopy
- 2x 1 W LED
- With two flexible gooseneck arms
- 2 years warranty



● LE.5208

TECHNICAL SPECIFICATIONS

- Economical dual LED intensity adjustable light source
- Equipped with two 44 cm flexible gooseneck arms, each equipped with a 1 W LED
- The light source version does not require bulb changing like halogen light sources and has a longer life span and low energy consumption
- The lightsource produces white light and there is no color temperature shifting when dimming
- Ideal for life and materials sciences, various industrial applications and for stereo microscopy

MODEL

Power	2x 1 W LED (equivalent to 2x 20 W halogen)
Luminous flux	18,000 lux max. (2 LED's)
Color temperature	6,500 K
Operating voltage	100-240 V AC / 12 V DC (50/60 Hz)
Dimensions	(h) 70 x (w) 120 x (d) 170 mm
Weight	1.6 kg
Product Number	LE.5208