

## Cold light source compact single 20 W

### HIGHLIGHTS

- Compact cold light source sources
- For Materials Science en Stereo microscopy
- Halogen 20 W version
- With  $\varnothing$  6 mm fibre light conductor
- 2 years warranty

### TECHNICAL SPECIFICATIONS

- Very compact and robust cold light source
- Supplied with a 40 cm self-sustaining  $\varnothing$  6 mm glass fiber light conductor
- Equipped with a 20 W halogen lamp
- Light intensity is factory preset and can not be adjusted
- For 230 Volt operation

### ILLUMINATION

Equipped with a 20 W halogen lamp. Light intensity is factory preset and can not be adjusted

### MODEL

<b>Power</b>	20 W halogen
<b>Color temperature</b>	3,000 K
<b>Operating voltage</b>	110 V AC and 230 V AC versions (50/60 Hz)
<b>Dimensions</b>	(h) 80 x (w) 110 x (d) 165 mm
<b>Weight</b>	2.0 kg
<b>Product Number</b>	LE.5209



● LE.5209