

Magnifiers & Measuring Microscopes

HIGHLIGHTS

- High-quality magnifiers
- Excellent optical quality
- Measuring microscopes are built with a microscope objective for magnification from 20 to 100x times
- To observe all details of the object
- Used by biologists, jewelers, watchmakers and philatelists

MAGNIFIERS WITH HANDLES

- EPB-5020** Magnifier 4x, lens Ø 50 mm and handle Ø 8 mm. Fits on the EPB-5022
- EPB-5021** Magnifier 10x, aplanatic lens Ø 25 mm and handle Ø 8 mm. Fits on the EPB-5022
- EPB-5022** Heavy stand with adjustable cross clamp, base Ø 70 mm. Weight 450 grams (For EPB-5022 and EPB-5021)



FOLDING MAGNIFIERS

- EPB-5016** Aplanatic magnifier 6x, doublet. Lens Ø 20 mm. Field of view 25 mm
- EPB-5018** Aplanatic magnifier 15x, doublet. Lens Ø 18 mm. Field of view 14 mm
- EPB-5032** Achromatic magnifier 10x, triplet. Lens Ø 20 mm. Field of view 15 mm
- EPB-5033** Achromatic magnifier 20x, triplet. Lens Ø 20 mm. Field of view 8 mm
- EPB-5034** Achromatic magnifier 10x, triplet. Lens Ø 16 mm. Field of view 15 mm
- EPB-5036** Aplanatic magnifier 10x, doublet. Lens Ø 20 mm. Field of view 20 mm
- EPB-5037** Achromatic diamond magnifier 10x, 4 lenses Ø 21 mm. Field of view 21 mm



**FOLDING MAGNIFIERS
WHITE LED ILLUMINATION**

- EPB-5034-LED** Achromatic magnifier 10x triplet. Lens Ø 21 mm. White LED illumination
- EPB-5018-LED** Achromatic magnifier 15x triplet. Lens Ø 21mm. White LED illumination
- EPB-5033-LED** Achromatic magnifier 20x, triplet. Lens Ø 21 mm. White LED illumination
- EPB-5035-LED** Achromatic magnifier 30x triplet. Lens Ø 21 mm. White LED illumination



**FOLDING MAGNIFIERS
WHITE LED AND UV ILLUMINATION**

- EPB-5034-LUV** Achromatic magnifier 10x triplet. Lens Ø 21 mm. White LED and UV illumination
- EPB-5018-LUV** Achromatic magnifier 15x triplet. Lens Ø 21mm. White LED and UV illumination
- EPB-5033-LUV** Achromatic magnifier 20x, triplet. Lens Ø 21 mm. White LED and UV illumination
- EPB-5035-LUV** Achromatic magnifier 30x triplet. Lens Ø 21 mm. White LED and UV illumination



TABLE MAGNIFIERS

- EPB-5039** Table magnifier, 15x magnification.
Glass achromatic lens Ø 19 mm. Field of view 15 mm
- EPB-5040** Table magnifier, 8x magnification. Plastic lens Ø 23 mm.
Field of view 30 mm. 10mm/100 reticle
- EPB-5041** Table magnifier, 10x magnification.
Glass lens Ø 26 mm. Field of view 25 mm
- EPB-5043** Table magnifier, 4x magnification.
Plastic lens Ø 57 mm. Field of view 55 mm




MEASURING MAGNIFIERS WITH RETICLE

- EPB-5044** Measuring magnifier 7x. Eyepiece Ø 18 mm. Field of view 24 mm.
Measuring scale 10-0-10 mm divided into 0.1 mm and crosshairs
- EPB-5045** Measuring magnifier 7x. Eyepiece Ø 18 mm. Field of view 24 mm.
Measuring scale 10-0-10 mm divided into 0.1 mm and crosshairs.
With 5 interchangeable scales for angle, thickness and linear measurements (see image next page)
- EPB-5046** Measuring magnifier 10x. Eyepiece Ø 20 mm. Field of view 35 mm.
Measuring scale 15-0-15 mm with 0.1, 0.5 and 1.0 mm graduations
- EPB-5047** Measuring magnifier 10x with LED. Eyepiece Ø 25 mm.
Field of view 30 mm. Measuring scale 15-0-15 mm divided into 0.1, 0.5 and 1.0 mm graduations. White LED illumination
- EPB-5048** Measuring magnifier 4x. Eyepiece Ø 30 mm. Field of view 65 mm.
Measuring scale 25-0-25 mm with 0.1, 0.5 and 1.0 mm graduations





● EPB-5045



MISCELLANEOUS MAGNIFIERS

EPB-5030 Dual sliding magnifier with 3x and 6x magnification.
Lens aperture 23 x 23 mm. Plastic housing.
Fields of view 50 x 50 mm and 25 x 25 mm

EPB-5049 Linen tester folding magnifier with 6x aplanatic lens.
A 25 x 25 mm base with millimeter scale.
Metal casing lacquered in black. Supplied with plastic case

MEASURING MICROSCOPES

EPB-5008 Measuring microscope 20x with LED illumination.
Field of view 7.2 mm

EPB-5010 Measuring microscope 40x with LED illumination.
Field of view 3.8 mm

EPB-5012 Measuring microscope 60x with LED illumination.
Field of view 1.45 mm

EPB-5014 Measuring microscope 100x with LED illumination.
Field of view 0.7 mm



● EPB-5030

● EPB-5049



EPB-5008 ●



www.euromex.com