

Miscellaneous accessories

ABOUT

- Stain through
- Fixation
- Cleaning products

PRODUCTS

- EPB-5200** Glass Staining Dish with glass cover, holds 10 microscope slides
- EPB-5210** Fine pointed brush
- EPB-5245** Lens cleaning paper, 100 sheets per pack
- EPB-5250** Solid paraffin 200 gram
- EPB-5255** Immersion oil, 25 ml. Refraction index $n = 1.515$
- EPB-5265** Entellan, quick drying Canada balsam. In bottle (25 ml)
- EPB-5274** Iso propyl alcohol 99%, 200 ml
- EPB-5277** Formalin 40%, fixing agent. In bottle (200 ml)
- EPB-5275** Cleaning kit: lens fluid, lint free lens tissue paper, brush, air blower, cotton swabs
- EPB-5276** Microscope maintenance and servicing kit, 16pcs: cleaning brush, 6 pcs screwdriver set, air blower, 3 pcs Allen key, 1.5, 2, 2.5 mm, lens cleaning fluid 20 ml, cleaning cloth 140 x 140 mm, 100 pcs Lens tissue sheets, tube of maintenance grease, 10 ml bottle of oil, packed in a nice toolbox



● EPB-5200



● EPB-5245



● EPB-5210



● EPB-5250



● EPB-5255



● EPB-5265



● EPB-5275



● EPB-5276



● EPB-5274



● EPB-5277

