

Preparation accessoires

HIGHLIGHTS

- Euromex offers the right tools for preparing samples
- Separate instruments and complete dissection kits
- Superior quality
- All our instruments are made from stainless steel

MODELS

DISSECTION INSTRUMENTS

- PB.5050** Dissecting needle, straight tip with plastic handle
PB.5054 Dissecting needle, lancet-shaped with plastic handle
PB.5056 Blunt-pointed probe, Ø2 mm, length 140 mm
PB.5058 Section lifter with wooden handle
- PB.5060** Serrated thumb dressing forceps, length 130 mm
PB.5062 Dressing forceps with blunt tips, length 130 mm
PB.5064 Dressing forceps with blunt tips, length 200 mm
PB.5067 Dissecting forceps with sharp tips, length 130 mm
PB.5069 Dissecting forceps with finely curved sharp tips, length 130 mm
PB.5077 Cover glass forceps, length 110 mm Scissors and blades



- PB.5080** Dressing scissors with blunt/sharp tip, length 140 mm
PB.5084 Dissecting scissors, sharp tips, length 110 mm
PB.5092 Scalpel with blunt tip, length 150 mm
PB.5094 Scalpel handle no. 4 (for blades PB.5095 and PB.5096)
PB.5095 Scalpel blades with sharp tip, five pieces per set
PB.5096 Scalpel blades with blunt tip, five pieces per set
PB.5103 Strop, double sided for stropping of scalpels, microtome knives. Includes stropping paste
PB.5105 Razor plano/concave, for sectioning. With plastic handle
PB.5106 Microtome knife plano/concave. With fixed handle, length 170 mm



HANDHELD MICROTOMES

- MT.5500** Cylinder microtome. Round stage Ø 70 mm with flat black glass top. Height adjustable with micrometer division 0.01 mm. Complete with plano-concave knife PB.5106 and wooden cabinet
- MT.5501** Cylinder microtome. Identical to MT.5500 but without knife and wooden cabinet
- MT.5503** Handheld microtome. Box type with table clamp and razor blade holder. Thickness adjustment by 0.025 mm intervals. Complete with wooden cabinet



DISSECTION KITS

PB.5110 4-PC Dissection kit in imitation leather case

Kit with dissecting needle straight (PB.5050), dissecting forceps (PB.5067), dissecting scissors (PB.5084) and scalpel (PB.5092).

Supplied with imitation leather case

PB.5110 ●



PB.5111 6-PC Dissection kit in plastic case

Kit with dissecting needle straight (PB.5050), dissecting needle lancet-shaped (PB.5054), dissecting forceps (PB.5067), dissecting scissors (PB.5084), scalpel handle no. 4 (PB.5094) and blades (PB.5096). Supplied with plastic case

● PB.5111



PB.5112 9-PC Dissection kit in imitation leather case

Kit with two straight dissecting needles (PB.5050), lancet-shaped dissecting needle (PB.5054), dressing forceps (PB.5062), dissecting forceps (PB.5067), scissors (PB.5080), dissecting scissors (PB.5084), two scalpels (PB.5094), five scalpel blades with straight tip and five blades with blunt tip (PB.5095/PB.5096).

Supplied with imitation leather case

PB.5112 ●



PB.5114 9-PC Dissection kit in wooden cabinet

Kit with two straight dissecting needles (PB.5050), lancet-shaped dissecting needle (PB.5054), dressing forceps (PB.5062), dissecting forceps (PB.5067), scissors (PB.5080), dissecting scissors (PB.5084), two scalpels (PB.5094), five scalpel blades with sharp tip (PB.5095) and five blades with blunt tip (PB.5096).

Supplied with wooden cabinet

● PB.5114



PB.5115 8-PC Dissection kit for Koi fish in plastic case

Dissecting set for Koi fish. Consisting of plastic box with forceps with 1 x 2 points (PB.5060), forceps blunt points, 130 mm (PB.5062), scissors, blunt/ sharp points, 140 mm (PB.5080), microscope scissors, sharp points, 110 mm (PB.5084), blunt-pointed probe (PB.5056), scalpel handle no. 4 (PB.5094), five scalpel blades with straight tip and five blades with blunt tip (PB.5095/PB.5096)

KIT FOR KOI FISH

PB.5115 ●



PB.5125 10-PC Professional slide preparation kit in plastic deluxe case

Kit with three prepared slides, 50 microscope slides (PB.5150) and 100 cover glasses (PB.5165). Two microscope slides with cavity. Canada balsam, dissecting needle (PB.5050 and PB.5054). Dissecting forceps (PB.5067), dissecting scissors (PB.5084) and scalpel (PB.5092). Supplied with plastic deluxe case

● PB.5125



All our superior quality instruments are made from stainless steel

