

Novex Explorer

HIGHLIGHTS

- High value for money
- 20x magnification
- Portable
- 2 year warranty

The Novex Explorer is the entry level stereo microscope for hobbyists, children and lower education. This light weight, cost effective, microscope enables users to see any object from stamp to ant at 20x magnification in bright crisp image quality. Due to its small size and weight the Novex explorer is ideal to take into the field. The LED's produce little heat, making them safe for live specimen and fingers of users. The ideal product as a first serious step into microscopy

TECHNICAL SPECIFICATIONS

OCULARS

Pair of WF10x /16mm eyepieces

OBJECTIVE

- Stereo objective 2x (fixed)
- 20 x magnification

STAND

Ergonomic carrying handle

ILLUMINATION

- Incident LED illumination
- Runs on batteries, (2x AA), not included

PACKAGE CONTENT

Supplied in Styrofoam case with eye cups and dust cover



● Explorer

NOVEX EXPLORER
220 (h) x 100 (w) x 135 (d) | 0.5 kg

